



**SCIENTIFIC
LABORATORY
SUPPLIES**

Material description: ETHYL BENZOATE pure
CAS No: 93-89-0

Supplied as: Product code Description
CHE1970 ETHYL BENZOATE pure 100ml.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99			
Boiling range	°C	211	213		
Wt. per ml.	g	1.051	1.052		

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1097/MS/19890101