

Material description: FEHLING'S SOLUTION NO. 2 QUALITATIVE  
CAS No: Mixture

<i>Supplied as:</i>	Product code	Description
	CHE1984	FEHLING'S SOLUTION NO. 2 QUALITATIVE 1L.
	CHE1986	FEHLING'S SOLUTION NO. 2 QUALITATIVE 500ml.
	CHE1988	FEHLING'S SOLUTION NO. 2 QUALITATIVE 2.5L.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Functional test <i>by Fehling's test</i>					Red precipitate formed
Density @20°C	g/ml	1.248	1.268	1.258	

Revision date: 30/04/2020

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1120/MS/20200430