



Material description: GRAM'S IODINE SOLUTION

CAS No: Mixture

<i>Supplied as:</i>	Product code	Description
	CHE2086	GRAM'S IODINE SOLUTION 100ml.
	CHE2088	GRAM'S IODINE SOLUTION 250ml.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Strength	mol/L	0.01280	0.01295	0.01287	

Revision date: 09/08/2012

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 3093/MS/20120809