

Material description: HEXAN-1-OL pure
CAS No: 111-27-3

Supplied as: Product code Description
CHE2102 HEXAN-1-OL pure 500ml.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	98			
Boiling range	°C	156	157		
Wt. per ml.	g	0.813	0.815		

Revision date: 15/10/2020

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0350/MS/20201015