



Material description: HYAMINE 1622 0.004M
CAS No: Mixture

Supplied as: Product code Description
 CHE2108 HYAMINE 1622 0.004M 2.5L.
 CHE2110 HYAMINE 1622 0.004M 1L.

Specification

| Parameter | Units | Minimum value | Maximum value | Typical value | Notes |
|-----------|-------|---------------|---------------|---------------|-------|
| Strength | mol/L | 0.003992 | 0.004008 | 0.004000 | |

Revision date: 19/04/2018 Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2385/MS/20180419