



Material description: HYDROCHLORIC ACID 3.5M (3.5N)

CAS No: 7647-01-0

Supplied as: Product code Description
CHE2130 HYDROCHLORIC ACID 3.5M (3.5N) 1L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Strength	mol/L	3.4965	3.5035		

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2564/MS/19890101