

Material description: HYDROCHLORIC ACID 35% w/w pure
CAS No: 7647-01-0

<i>Supplied as:</i>	Product code	Description
	CHE2136	HYDROCHLORIC ACID 35% w/w pure 2.5L.
	CHE2138	HYDROCHLORIC ACID 35% w/w pure 1L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Strength	% w/w	35	38		
Wt. per ml.	g			1.18	
Free chlorine	%		0.0004		

Revision date: 09/05/2018

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0361/MS/20180509