



Material description: HYDROCHLORIC ACID 35% w/w (in safety coated glass)

CAS No: 7647-01-0

Supplied as: Product code Description
CHE2140 HYDROCHLORIC ACID 35% w/w (in safety coated glass) 2.5L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Strength	% w/w	35	38		
Wt. per ml.	g			1.18	
Free chlorine	%		0.0004		

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2397/MS/19890101