



Material description: HYDROCHLORIC ACID 5% w/v
CAS No: 7647-01-0

Supplied as: Product code Description
CHE2150 HYDROCHLORIC ACID 5% w/v 10L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Strength	mol/L	1.33	1.41	1.37	

Revision date: 15/12/2009

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 3185/MS/20091215