

Material description: HYDROCHLORIC ACID 5.0M (5N)

CAS No: 7647-01-0

<i>Supplied as:</i>	Product code	Description
	CHE2156	HYDROCHLORIC ACID 5.0M (5N) 1L.
	CHE2158	HYDROCHLORIC ACID 5.0M (5N) 2.5L.
	CHE2160	HYDROCHLORIC ACID 5.0M (5N) 5L.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Strength	mol/L	4.995	5.005		

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2521/MS/19890101