

Material description: HYDROCHLORIC ACID 1.0M (N/1)

CAS No: 7647-01-0

<i>Supplied as:</i>	Product code	Description
	CHE2164	HYDROCHLORIC ACID 1.0M (N/1) 2.5L.
	CHE2166	HYDROCHLORIC ACID 1.0M (N/1) 1L.
	CHE2168	HYDROCHLORIC ACID 1.0M (N/1) 5L.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Strength	mol/L	0.999	1.001		

Revision date: 02/03/2022

Issue date: 12/05/2022

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1735/MS/20220302