

Material description: `4.5' INDICATOR SOLUTION

CAS No: Mixture

Supplied as: Product code Description
CHE2222 `4.5' INDICATOR SOLUTION 500ml.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Colour at pH 3.5; Orange / Red					ph 3.5
Colour at pH 4.4; Orange / Grey					pH 4.4
Colour at pH 4.5; Grey					pH 4.5
Colour at pH 4.6; Blue/Grey					pH 4.6
Colour at pH 6.0; Deep blue					pH 6.0

Revision date: 26/02/2020

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2289/MS/20200226