



Material description: IODINE A.R.
CAS No: 7553-56-2

Supplied as: Product code Description
CHE2236 IODINE A.R. 100g.
CHE2238 IODINE A.R. 250g.

Specification

| Parameter | Units | Minimum value | Maximum value | Typical value | Notes |
|---------------|-------|---------------|---------------|---------------|-------|
| Purity | % | 99.8 | | | |
| Non-volatiles | % | | 0.02 | | |

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0377/MS/19890101