

Material description: LACTIC ACID pure
CAS No: 79-33-4

<i>Supplied as:</i>	Product code	Description
	CHE2282	LACTIC ACID pure 250ml.
	CHE2284	LACTIC ACID pure 500ml.
	CHE2286	LACTIC ACID pure 2.5L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Total acidity-lactic	% w/w	79.5	80.5	80.0	
Wt. per ml.	g	1.18	1.20	1.18	
Colour - Fresh solution	Hazen			50	

Revision date: 04/08/2016

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1185/MS/20160804