



Material description: LEAD METAL FOIL

CAS No: 7439-92-1

Supplied as: Product code Description
CHE2322 LEAD METAL FOIL 500g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity; as Lead	%	99			
Nominal thickness	mm			0.45	

Revision date: 10/07/2009

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1606/MS/20090710