



Material description: LITMUS SOLUTION

CAS No: Mixture

<i>Supplied as:</i>	Product code	Description
	CHE2370	LITMUS SOLUTION 100ml.
	CHE2372	LITMUS SOLUTION 500ml.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Colour change; Red to Blue	pH				pH 5.0-8.0

Revision date: 18/02/2020

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1213/MS/20200218