



**SCIENTIFIC
LABORATORY
SUPPLIES**

Material description: L-LYSINE MONOHYDROCHLORIDE

CAS No: 657-27-2

Supplied as: Product code

CHE2374

Description

L-LYSINE MONOHYDROCHLORIDE 25g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Assay	%	98.0	102.0		

Revision date: 25/05/2021

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 3298/MS/20210525