

Material description: 2-METHYLPROPAN-1-OL pure
CAS No: 78-83-1

Supplied as:

Product code	Description
CHE2400	2-METHYLPROPAN-1-OL pure 500ml.
CHE2402	2-METHYLPROPAN-1-OL pure 2.5L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	98			
Boiling range	°C	105	109		
Wt. per ml.	g	0.800	0.803		
Water	% w/w		0.2		
Non-volatiles	%		0.003		

Revision date: 18/10/1993 Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1261/MS/19931018