

Material description: 2-METHYLPROPAN-1-OL pure  
CAS No: 78-83-1

*Supplied as:* Product code Description  
CHE2400 2-METHYLPROPAN-1-OL pure 500ml.  
CHE2402 2-METHYLPROPAN-1-OL pure 2.5L.

## Specification

| Parameter     | Units | Minimum value | Maximum value | Typical value | Notes |
|---------------|-------|---------------|---------------|---------------|-------|
| Purity        | %     | 98            |               |               |       |
| Boiling range | °C    | 105           | 109           |               |       |
| Wt. per ml.   | g     | 0.800         | 0.803         |               |       |
| Water         | % w/w |               | 0.2           |               |       |
| Non-volatiles | %     |               | 0.003         |               |       |

Revision date: 18/10/1993 Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1261/MS/19931018