

Material description: METHANOL A.R.

CAS No: 67-56-1

<i>Supplied as:</i>	Product code	Description
	CHE2410	METHANOL A.R. 1L.
	CHE2528	METHANOL A.R. 2.5L.
	CHE2530	METHANOL A.R. 25L.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.85			
Water	% w/w		0.1		
Acidity; acetic	%		0.003		
Iron	mg/Kg		0.1		
Chloride; as Cl	mg/Kg		0.5		
Sulphur	mg/Kg		0.5		
Carbonilic comp.	mg/Kg		30		
Non-volatiles	mg/L		8		
Colour	Hazen		5		

Revision date: 30/10/2009

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0428/MS/20091030