

Material description: 4-METHYLPENTAN-2-ONE pure
CAS No: 108-10-1

<i>Supplied as:</i>	Product code	Description
	CHE2412	4-METHYLPENTAN-2-ONE pure 500ml.
	CHE2414	4-METHYLPENTAN-2-ONE pure 2.5L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99			
Boiling range	°C	114	117		
Wt. per ml.	g	0.799	0.802		
Water	% w/w		0.1		

Revision date: 16/09/1994

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0438/MS/19940916