

Material description: MAGNESIUM CARBONATE BASIC HEAVY

CAS No: 39409-82-0

<i>Supplied as:</i>	Product code	Description
	CHE2418	MAGNESIUM CARBONATE BASIC HEAVY 250g.
	CHE2420	MAGNESIUM CARBONATE BASIC HEAVY 500g.
	CHE5201	MAGNESIUM CARBONATE BASIC HEAVY 25Kg.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity ; as MgO	%	40	45		
Arsenic	%		0.002		
Iron	%		0.04		
Lead	%		0.0001		
Water soluble matter	%		1		

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1609/MS/19890101