

Material description: METHANOL- for Karl Fisher & non aqueous titrations
CAS No: 67-56-1

<i>Supplied as:</i>	Product code	Description
	CHE2524	METHANOL- for Karl Fisher & non aqueous titrations 500ml.
	CHE2526	METHANOL- for Karl Fisher & non aqueous titrations 2.5L.
	CHE4006	METHANOL- for Karl Fisher & non aqueous titrations 1L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Wt. per ml.	g	0.790	0.792		
Water	% w/w		0.01		

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1908/MS/19890101