

Material description: METHANOL HPLC
CAS No: 67-56-1

Supplied as: Product code Description
CHE2532 METHANOL HPLC 2.5L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.8			
Acidity; acetic	meq/g		0.0002		
Water	% w/w		0.05		
Non-volatiles	ppm		5.0		
Absorbance @ 210nm	AU		0.200		
Absorbance @ 220nm	AU		0.100		
Absorbance @ 230nm	AU		0.050		
Absorbance @ 240nm	AU		0.020		
Absorbance @ 250nm	AU		0.005		
Absorbance @ 260nm	AU		0.005		

Revision date: 29/04/2021 Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2128/MS/20210429