

Material description: METHANOL pure  
CAS No: 67-56-1

<i>Supplied as:</i>	Product code	Description
	CHE2534	METHANOL pure 2.5L.
	CHE2536	METHANOL pure 1L.
	CHE5010	METHANOL pure 25L.
	CHE5026	METHANOL pure 50L.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99			
Wt. per ml.	g	0.791	0.793		

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0427/MS/19890101