

Material description: METHYL BENZOATE pure
CAS No: 93-58-3

Supplied as: Product code Description
CHE2548 METHYL BENZOATE pure 500ml.
CHE2550 METHYL BENZOATE pure 100ml.

Specification

| Parameter | Units | Minimum value | Maximum value | Typical value | Notes |
|-------------------------|-------|---------------|---------------|---------------|-------|
| Purity | % | 99 | | | |
| Boiling point; @10mm Hg | °C | | | 84 | |
| Wt. per ml. | g | 1.086 | 1.088 | | |

Revision date: 28/11/1991

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2352/MS/19911128