



Material description: METHYL ORANGE (C.I. 13025)

CAS No: 547-58-0

<i>Supplied as:</i>	Product code	Description
	CHE2558	METHYL ORANGE (C.I. 13025) 25g.
	CHE2560	METHYL ORANGE (C.I. 13025) 50g.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Orange red-orange yellow	pH	2.9	4.6		

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1844/MS/19890101