



**SCIENTIFIC  
LABORATORY  
SUPPLIES**

Material description: METHYL ORANGE SOLUTION

CAS No: Mixture

<i>Supplied as:</i>	Product code	Description
	CHE2562	METHYL ORANGE SOLUTION 500ml.
	CHE2564	METHYL ORANGE SOLUTION 100ml.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Colour change; Orange red to Orange yellow					pH 2.9-4.6

Revision date: 09/05/2018

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1258/MS/20180509