



Material description: METHYL ORANGE SOLUTION (screened)

CAS No: Mixture

<i>Supplied as:</i>	Product code	Description
	CHE2566	METHYL ORANGE SOLUTION (screened) 500ml.
	CHE2568	METHYL ORANGE SOLUTION (screened) 100ml.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Colour change; Violet to Green					pH 2.9-4.6

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1845/MS/19890101