



**SCIENTIFIC
LABORATORY
SUPPLIES**

Material description: METHYL RED SOLUTION

CAS No: Mixture

<i>Supplied as:</i>	Product code	Description
	CHE2572	METHYL RED SOLUTION 100ml.
	CHE2574	METHYL RED SOLUTION 500ml.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Colour change; Red to Yellow					pH 4.2-6.3

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1263/MS/19890101