

Material description: METHYLENE BLUE 1% w/v
CAS No: Mixture

<i>Supplied as:</i>	Product code	Description
	CHE2586	METHYLENE BLUE 1% w/v 100ml.
	CHE2588	METHYLENE BLUE 1% w/v 250ml.
	CHE3958	METHYLENE BLUE 500ml.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Absorbance @ 662nm (0.1%)	AU	1.850	1.950		

Revision date: 15/08/2019 Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 3153/MS/20190815