



Material description: MOLECULAR SIEVE TYPE 3A (1/16 inch pellets)

CAS No:

<i>Supplied as:</i>	Product code	Description
	CHE2596	MOLECULAR SIEVE TYPE 3A (1/16 inch pellets) 500g.
	CHE2598	MOLECULAR SIEVE TYPE 3A (1/16 inch pellets) 1Kg.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Typical pore diameter	nm			0.3	

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0445/MS/19890101