

Material description: N,N-DIMETHYLFORMAMIDE pure
CAS No: 68-12-2

<i>Supplied as:</i>	Product code	Description
	CHE2602	N,N-DIMETHYLFORMAMIDE pure 2.5L.
	CHE2604	N,N-DIMETHYLFORMAMIDE pure 500ml.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99			
Water	% w/w		0.03		
Colour	Hazen		10		
Acidity; acetic	%		0.001		
Base; dimethylamine	%		0.001		
Methanol	%		0.005		

Revision date: 01/11/2007

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0278/MS/20071101