

Material description: n-BUTYL ACETATE pure
CAS No: 123-86-4

<i>Supplied as:</i>	Product code	Description
	CHE2612	n-BUTYL ACETATE pure 500ml.
	CHE2614	n-BUTYL ACETATE pure 2.5L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	98			
Boiling range	°C	124	127		
Wt. per ml.	g	0.878	0.882		
Acidity; acetic	%		0.05		
Water	% w/w		0.1		

Revision date: 01/08/1997 Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0128/MS/19970801