

Material description: n-HEPTANE HPLC
CAS No: 142-82-5

Supplied as: Product code Description
CHE2626 n-HEPTANE HPLC 2.5L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99			
Water	% w/w		0.010		
Absorbance @ 200nm	AU		0.250		
Absorbance @ 210nm	AU		0.100		
Absorbance @ 220nm	AU		0.050		
Absorbance @ 230nm	AU		0.020		
Absorbance @ 240nm	AU		0.010		
Absorbance @ 250nm	AU		0.005		
Residue on evaporation	ppm		5.00		

Revision date: 29/06/2020

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2933/MS/20200629