

Material description: n-HEXANE (fraction from petroleum)

CAS No: 110-54-3

<i>Supplied as:</i>	Product code	Description
	CHE2632	n-HEXANE (fraction from petroleum) 2.5L.
	CHE2634	n-HEXANE (fraction from petroleum) 500ml.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	95			
Boiling range	°C	67	69		
Wt. per ml.	g			0.66	

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2620/MS/19890101