

Material description: n-HEXANE (fraction from petroleum) HPLC
CAS No: 110-54-3

Supplied as: Product code Description
CHE2636 n-HEXANE (fraction from petroleum) HPLC 2.5L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Assay	%	99.00			
Water	ppm		30		
Wt. per ml.	g			0.665	
Transmission; 250nm	%	99			
Transmission; 240nm	%	98			
Transmission; 230nm	%	95			
Transmission; 220nm	%	87			
Transmission; 210nm	%	63			

Revision date: 14/04/2021

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2283/MS/20210414