



**SCIENTIFIC
LABORATORY
SUPPLIES**

Material description: n-HEXANE pure
CAS No: 110-54-3

Supplied as:	Product code	Description
	CHE2638	n-HEXANE pure 500ml.
	CHE2640	n-HEXANE pure 100ml.
	CHE2642	n-HEXANE pure 2.5L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	95			
Boiling range	°C	68	69		
Density @15°C	g/ml	0.663	0.669		

Revision date: 02/03/2021

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2167/MS/20210302