

Material description: NICKEL (II) SULPHATE 6H₂O pure
CAS No: 10101-97-0

<i>Supplied as:</i>	Product code	Description
	CHE2660	NICKEL (II) SULPHATE 6H ₂ O pure 250g.
	CHE2662	NICKEL (II) SULPHATE 6H ₂ O pure 500g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity; as Nickel	%			22	
Cobalt	%		0.02		
Iron	%		0.005		
Insoluble matter	%		0.05		

Revision date: 21/11/2002

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0465/MS/20021121