

Material description: NITRIC ACID 69% w/w A.R.

CAS No: 7697-37-2

Supplied as: Product code Description
CHE2670 NITRIC ACID 69% w/w A.R. 2.5L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Strength	% w/w	68	70		
Wt. per ml.	g			1.42	
Sulphate	ppm		2		
Arsenic	ppm		0.010		
Iron	ppm		0.20		
Lead	ppm		0.050		
Non-volatiles	ppm		10		

Revision date: 09/08/2021

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0469/MS/20210809