



Material description: NITRIC ACID 0.1M (N/10)

CAS No: 7697-37-2

Supplied as: Product code Description  
CHE2684 NITRIC ACID 0.1M (N/10) 2.5L.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Strength	mol/L	0.0999	0.1001		

Revision date: 09/05/2018

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1764/MS/20180509