

Material description: n-OCTANE pure  
CAS No: 111-65-9

<i>Supplied as:</i>	Product code	Description
	CHE2686	n-OCTANE pure 100ml.
	CHE2688	n-OCTANE pure 500ml.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99			
Boiling range	°C	124	126		
Wt. per ml.	g	0.701	0.703		

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2160/MS/19890101