

Material description: n-PENTANE pure  
CAS No: 109-66-0

*Supplied as:* Product code Description  
CHE2694 n-PENTANE pure 2.5L.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	98			
Boiling range	°C	35	36		
Wt. per ml.	g	0.625	0.627		

Revision date: 28/01/2021 Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2164/MS/20210128