



**SCIENTIFIC
LABORATORY
SUPPLIES**

Material description: PARAFFIN LIQUID (heavy)

CAS No: 8042-47-5

Supplied as: Product code Description
CHE2740 PARAFFIN LIQUID (heavy) 2.5L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Wt. per ml.	g			0.87	
Kinematic viscosity @40°C	cSt	63	83		

Revision date: 01/08/1990

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2387/MS/19900801