

Material description: PARAFFIN LIQUID (medium)

CAS No: 8042-47-5

<i>Supplied as:</i>	Product code	Description
	CHE2744	PARAFFIN LIQUID (medium) 2.5L.
	CHE2746	PARAFFIN LIQUID (medium) 25L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Wt. per ml.	g	0.81	0.82		
Kinematic viscosity @40°C	cSt	16	20		

Revision date: 01/08/1990

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1288/MS/19900801