

Material description: PENTAN-1-OL (n-Amyl alcohol)

CAS No: 71-41-0

*Supplied as:* Product code Description  
CHE2758 PENTAN-1-OL (n-Amyl alcohol) 250ml.

## Specification

| Parameter     | Units | Minimum value | Maximum value | Typical value | Notes |
|---------------|-------|---------------|---------------|---------------|-------|
| Purity        | %     | 98            |               |               |       |
| Boiling range | °C    | 136           | 138           |               |       |
| Wt. per ml.   | g     | 0.814         | 0.815         |               |       |

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1853/MS/19890101