

Material description: PERCHLORIC ACID 60% w/w pure
CAS No: 7601-90-3

Supplied as: Product code Description
CHE2768 PERCHLORIC ACID 60% w/w pure 500ml.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Strength	% w/w	60	62		
Wt. per ml.	g			1.54	
Chloride	%		0.005		
Iron	%		0.0001		
Lead	%		0.0001		
Non-volatiles	%		0.005		

Revision date: 09/05/1990 Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0506/MS/19900509