



Material description: PERIODIC ACID pure  
CAS No: 10450-60-9

Supplied as:	Product code	Description
	CHE2776	PERIODIC ACID pure 50g.
	CHE2778	PERIODIC ACID pure 100g.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	98			
Melting range	°C	124	129		

Revision date: 29/10/1996 Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1772/MS/19961029