

Material description: PETROLEUM SPIRIT 100/120 pure

CAS No: 8032-32-4

| <i>Supplied as:</i> | Product code | Description                          |
|---------------------|--------------|--------------------------------------|
|                     | CHE2782      | PETROLEUM SPIRIT 100/120 pure 2.5L.  |
|                     | CHE2784      | PETROLEUM SPIRIT 100/120 pure 1L.    |
|                     | CHE2786      | PETROLEUM SPIRIT 100/120 pure 500ml. |
|                     | CHE5060      | PETROLEUM SPIRIT 100/120 pure 25L.   |

## Specification

| Parameter          | Units | Minimum value | Maximum value | Typical value | Notes |
|--------------------|-------|---------------|---------------|---------------|-------|
| Boiling range; 90% | °C    | 95            | 125           |               |       |
| Wt. per ml.        | g     | 0.700         | 0.745         |               |       |

Revision date: 11/07/1990

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1305/MS/19900711