



**SCIENTIFIC
LABORATORY
SUPPLIES**

Material description: PETROLEUM SPIRIT 40/60 pure

CAS No: 8032-32-4

Supplied as:	Product code	Description
	CHE2796	PETROLEUM SPIRIT 40/60 pure 2.5L.
	CHE2798	PETROLEUM SPIRIT 40/60 pure 500ml.
	CHE2800	PETROLEUM SPIRIT 40/60 pure 25L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Boiling range; 90%	°C	40	60		
Wt. per ml.	g	0.655	0.670		

Revision date: 11/07/1990

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1302/MS/19900711